TREATMENT LED MASK FOR THE FACE

USER MANUAL

PROFESSIONAL TREATMENT AT HOME

DEAR CUSTOMER,

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EN

PALSAR7 - ABOUT THE BRAND

PALSAR7 IS A BRAND SPECIALIZING IN INNOVATIVE PRODUCTS FOR PROFESSIONAL SKIN CARE FROM THE COMFORT OF YOUR HOME. IT OFFERS MODERN TECHNOLOGIES THAT SUPPORT SKIN HEALTH AND BEAUTY WITHOUT THE NEED TO VISIT A BEAUTY SALON. THE BRAND EMPHASIZES QUALITY, SAFETY AND EASE OF USE OF ITS PRODUCTS. THANKS TO THIS, YOU CAN INDULGE IN PROFESSIONAL CARE ANYTIME, ANYWHERE.

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TREATMENT LED MASK FOR THE FACE

PHOTONE LIGHT THERAPY

Why go for expensive chemical peels and skin cleansers, laser treatments and hyaluronic acid wrinkle fillers when you can treat your skin from the comfort of your home?

The PALSAR7 treatment LED mask is an innovative way to get rid of acne, wrinkles, tired skin and other imperfections. This is thanks to photon light therapy, which is one of the most effective tools in the field of dermatology in the fight against chronic skin defects. Unlike laser treatments, photon treatment with LED lights does not cause any pain or side effects that require time to heal.

HOW PHOTONE THERAPY WORKS

Light can be simply described as a visible part of electromagnetic radiation of different frequencies. It reaches different wavelengths in the range from 390 to 760nm. In this therapy, light penetrates the body as a stream of energy in the form of photons. Different wavelengths of light have different positive effects on skin care.

LIGHT WAVELENGTHS FOR THE LED TREATMENT MASK:

Red light: **630 nm**

Yellow light: **590 nm**

Green light: **520 nm**

Blue light: **470 nm**

Using this mask is very simple, you just need to put it on, fasten it to your face with an elastic band and adjust the intensity and duration of the selected treatment with the controller. There are soft gray rubber pads around the eye holes to prevent unpleasant pressure around the eyes. The treatment LED face mask is equipped with 150 LED diodes that produce photons with different wavelengths. Each diode has a defined power and luminous flux value. Therefore, it is important to distinguish the individual colors of the diodes. One color treats acne, another treats wrinkles, scars, redness, swelling - it depends on your preference.

The LEDs on the mask do not need to be changed.

EFFECTS OF TREATING LED MASKS:

- Reduces wrinkles and prevents them from occurring
- Minimizes pores
- Fights acne
- Helps fight facial grease
- Improves skin elasticity
- Unifies skin tone
- Suppresses melanin production
- Improves blood circulation in subcutaneous tissue
- Increases the efficiency of absorption of active ingredients of cosmetics

MAIN FUNCTION:

- Natural light waves improve the appearance of the skin
- Each light color is effective for a specific skin treatment
- Different mode options
- Easy operation
- Easily portable

MASK DESCRIPTION:


THE MASK IS SUITABLE FOR THOSE WHO:

- Stay up late and their skin is tired.
- Sit in front of the computer screen or mobile screen for long periods of time.
- Do not follow a healthy diet that improves skin condition
- Have problematic skin: spots, acne, clogged pores, dull, yellowish, rough or sagging skin, etc.

INSTRUCTIONS FOR USE:

- **Before using the mask, test all functions of the controller without the mask on your face. Try to avoid direct eye contact with the light coming from the mask.**
- Connect the mask to the controller using the USB cable. Connect the controller to the power cable.
- **The power cord must be plugged in when using the mask.**
- Adjust the length of the elastic band to your face.
- Before putting on the mask, thoroughly clean the skin of the face and neck and then gently dry with a towel.
 - Press the **"power switch"** button on the controller. If the mask is on, the green **"power"** light on the top left will light up.
 - Then set the light intensity using the **"energy + and -"** buttons to your desired value.
 - Press the **"on/off"** button and the mask will light up. If the mask doesn't light up, you don't have the light intensity set. In this case, return to Step 2, otherwise continue with Step 4.
 - Use the **"color"** button to select the color of the treatment you want to undergo. Try not to look directly into the lights of the mask when choosing a color.
 - Use the **"time"** button to set the required treatment time, the recommended time is 10-15 minutes.
 - Put the mask on your face and start the whole skin treatment process with photon therapy.

CONTROLLER DESCRIPTION:

The maximum time of using the mask per day is 45 minutes. If you want to try several treatments in one day, wait at least 10 minutes between them. The recommended duration of one treatment is 10-15 minutes. Ideally, use the mask 2-3 times a week. If you suffer from sensitive skin, eczema, or allergies have a cold, consult your doctor before using the mask.

IMAGES WITH DESCRIPTION ARE FOR ILLUSTRATION PURPOSES ONLY.
MADE IN CHINA.

EFFECTS OF INDIVIDUAL LIGHTS:

1) 630NM-LED RED LIGHT: ANTI-WRINKLE

- Whitening effects
- For soft, wrinkle-free skin
- Smoothes fine wrinkles, shrinks pores
- Improves skin condition
- Supports the natural production of collagen
- Increases the absorption efficiency of active ingredients in cosmetic products
- Improves blood circulation in the subcutaneous tissue
- Prevents the formation of pigment

2) 520NM-LED GREEN LIGHT: ANTI-FATIGUE

- Improves skin condition
- Balances the color tone of the skin
- Reduces sebum production
- Relieves mental stress
- Reduces the occurrence of edema (swelling)
- Contributes to the elimination of epidermal pigment spots
- Slows down the production of melanin
- Has an anti-inflammatory effect (especially against the formation of acne skin)
- Suitable for all skin types

3) 470NM-LED BLUE LIGHT: AGAINST ACNE

- Effectively reduces the appearance and cures acne
- Improves the condition of the skin without leaving scars
- Activates porphyrins (they start a chemical reaction that destroys harmful bacteria)
- Has an antibacterial and anti-inflammatory effect

4) 590NM-LED YELLOW LIGHT: AGAINST REDNESS, TO INCREASE IMMUNITY

- Reduces skin roughness and redness
- Relieves skin diseases
- In addition, it relieves heat and helps to increase immunity
- Supports the absorption of cosmetics
- Contributes to the improvement of cellular metabolism
- Improves detoxification properties by supporting lymphatic activity
- Especially effective after intensive mechanical peeling
- Helps to slow down the aging of the skin

5) LED PURPLE LIGHT: AGAINST ACNE SCARS

- Combination of red and blue light = combination of two treatment therapies
- Particularly effective in treating acne scars

6) LED GREEN-BLUE LIGHT: METABOLISM SUPPORT

- Gradually increases cellular energy
- Supports metabolism

7) LED WHITE LIGHT: ANTI-AGING

- Treats age spots and sagging skin
- Reduces the appearance of fine wrinkles.

WARNING:

- Not suitable for children under 14 years of age.
- Forbidden to use: pregnant women, infants, children, people with colds, angina or flu, epileptics, people allergic to light, patients with thyroid disease.
- Forbidden to use by people who felt sick after using the mask - visit your doctor and consult with them about the use of the mask.
- When using the mask, do not watch TV, surf the Internet, read a book, etc. - there is a risk of long-term damage to the eyes.
- Avoid prolonged direct eye contact with the light source.
- If you intend to use cosmetics at the same time as the mask, choose high-quality products with a natural composition.
- After opening the box, the mask may have a slight smell, this is because it is made of medical silicone, the smell will disappear soon.
- Use the mask in bed or on a lounger, always keep it in the box to prevent damage.
- Do not use water to clean the mask, do not immerse the mask in water - use a damp cloth to clean.
- Do not use any chemical products to clean the mask.
- Do not expose the mask to direct sunlight.
- Use only as directed, be cautious.

ACCESSORIES:

1x LED face mask, 1x remote control, 1x USB cable, 1x charging adapter, 1x user manual

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Wear and tear of the product and its components caused by their common use (discharged batteries, loss of the original color of the product, decreased capacity of accumulator, wear and tear of mechanical or displaying components, wear and tear of supply cables, etc.) shall be subject neither to defect warranty nor quality warranty in the moment of take-over.

Furthermore, complaints about product parameters, not included in the operating instructions or in any commercial-technical documentation of Seller, related to the product, shall not be subject to quality warranty in the moment of take-over.

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