
LUCENT

PULSE

Water Flosser

MANUAL AND WARRANTY

PULSE

User Manual

INSTRUCTIONS OF USE	ENGLISH	4-5
HANDBUCH	DEUTSCH	6-7
MODE D'EMPLOI	FRANÇAIS	8-9
MANUAL DE USUARIO	ESPAÑOL	10-11
MANUALE D'USO	ITALIANO	12-13
GEBRUIKSAANWIJZING	NEDERLANDS	14-15
BRUKSANVISNING	SVENSKA	16-17
BRUGSANVISNING	NORSK	18-19
BRUGSANVISNING	DANSK	20-21
KÄYTTÖOHJEET	SUOMALAINEN	22-23
INSTRUKCJA UŻYCIA	POLSKI	24-25
NAVODILA ZA UPORABO	SLOVENŠČINA	26-27
NÁVOD K POUŽITÍ	ČEŠTINA	28-29
NÁVOD NA POUŽITIE	SLOVENČINA	30-31

THANK YOU FOR BUYING THE PULSE WATER FLOSSER

Meet Pulse — the water flosser designed to make deep cleaning feel effortless.

With targeted water pressure, three intuitive cleaning modes, and a compact design built for real bathrooms (and real travel bags), Pulse helps clean the places your toothbrush can't quite reach.

Use Pulse daily for healthier gums, a fresher mouth, and a dental routine that feels just a little smarter.

PRODUCT GUIDE AND SPECS

- 01** Nozzle
- 02** Water Inlet
- 03** Nozzle Button
- 04** Power Button
- 05** Mode Indicator
- 06** Mode Button
- 07** Water Tank

HOW TO USE PULSE

1. Fill the Water Tank

Open the water inlet (02) and fill the water tank (07) with room-temperature water. Do not use hot water. Close the inlet securely before use.

2. Attach the Nozzle

Insert the nozzle (01) into the top of the device until it clicks into place. To remove the nozzle, press the nozzle button (03) and pull the nozzle upward.

3. Turn On / Off

Press the power button (04) once to turn the device on or off.
4. Select Cleaning Mode

Press the mode button (06) to cycle through the three modes. The mode indicator lights (05) will show the selected mode.

5. Floss Your Teeth

Place the nozzle inside your mouth before turning the device on to avoid splashing. Lean slightly over the sink and guide the water stream along the gumline and between teeth. Move slowly from tooth to tooth for best results.

Avoid pressing the nozzle directly against teeth or gums.

6. After Use

Turn the device off and empty any remaining water from the tank. Allow the device and nozzle to dry before storing.

MODES

Soft

Gentle water pressure designed for sensitive teeth and first-time users.

Normal

Balanced water pressure for everyday cleaning and plaque removal.

Pulse

Pulsating water pressure designed for deeper cleaning and gum massage.

NOZZLES

Standard nozzles

For daily cleaning between teeth. Helps remove plaque and food debris.

Orthodontic nozzle

Designed to clean around braces, wires and other hard to reach areas.

Periodontal nozzle

Designed for gentle cleaning along the gumline and in periodontal pockets.

CLEANING AND MAINTENANCE

Cleaning the Water Tank

1. Remove the water tank (07) from the device.
2. Wash using mild soap and clean water.
3. Dry with a soft cloth before reattaching.

Regular cleaning helps maintain hygiene and performance.

Cleaning the Nozzle

Remove the nozzle (01) and rinse with clean water after each use. For optimal hygiene, replace the nozzle periodically.

PRECAUTIONS

- Do not operate the device without reading the instructions.
- Do not use hot water in the tank.
- Do not use the device if the charging cable is damaged.
- Do not charge the device in wet or damp areas.
- Do not store the device near heat sources or high humidity.
- Do not use the device for purposes other than oral cleaning.
- Do not press the nozzle directly against teeth or gums.
- Children should only use the device under adult supervision.
- Do not use the device without a nozzle attached.
- If you have had oral surgery or severe periodontal disease within the last two months, consult your dental professional before use.

TECHNICAL SPECIFICATIONS

Water tank capacity:	250 ml
Pressure range:	40–100 PSI
Pulse frequency:	up to 1,600 pulses per minute
Modes:	Soft / Normal / Pulse
Battery:	Rechargeable Li-ion battery
Battery capacity:	2500 mAh
Charging:	USB Type C
Input voltage:	5V and minimum 1A
Charging time:	approx. 4 hours
Battery life:	up to 30 days (based on typical use)
Waterproof rating:	IPX7

WARRANTY

Our products are designed and manufactured to high standards to ensure reliable performance. If a manufacturing defect occurs, we provide a two-year warranty from the date of purchase, covering defects in materials and workmanship. To claim warranty service, please present the original purchase receipt to the retailer or contact our support team. The warranty does not cover misuse, normal wear and tear, unauthorized repairs, or damage caused by external events.

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VÝROBCA**Be Lucent AB**