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The PALSAR7 Galvanic iron for the neck and décolleté is an innovative tool thanks to which you can clean and nourish your skin better than ever. The product cares for the skin with a thorough high-frequency circuit inside the iron, which makes the result very effective.

Before using the iron, always clean your neck and décolleté thoroughly with quality products and rinse it with warm water to achieve maximum effects.

If you do not clean your skin and the iron well, you can damage your skin. Please note that improper use of the iron and insufficient cleaning can lead to clogged pores, acne, rash, sagging skin, black spots and other skin problems.

What is the iron used for?

The iron has functional modes that help nourish and regenerate the skin, lift the skin, or smooth out wrinkles.

INSTRUCTIONS FOR USE

Charge the device for at least 2 hours before first use.

1. Before each use, clean the iron with a piece of clean cloth or a cotton swab.
2. Use only on clean skin. Works better with neck cream.
3. Press and hold the bottom button  for 2 seconds to turn the device on or off, then press the button once to switch the function.
4. After selecting the mode, the intensity can be adjusted using the top button. The selected mode and intensity is displayed on the screen.
5. During use, it is necessary to hold both metal strips on the sides of the handle.

POPIS PRODUKTU

TREATMENT MODES

Iron can be used daily for upto 5-10 minutes. The choice of mode is entirely up to you, you can choose your preferred one or switch between multiple modes within the time limit.

Skin rejuvenation

Skin rejuvenation mode emits blue light. It supports protein and collagen synthesis, which helps maintain skin firmness and reduces skin unevenness. Smoothes pores, helps sensitive skin and is great for oily skin.

Facial lifting

Skin lifting mode uses red light. Accelerates blood circulation and promotes collagen and fibroblast formation to smooth wrinkles and lines. Firms the skin. Lightens skin pigment and helps remove scars.

Wrinkle improver

Purple light is used in the wrinkle improvement mode. It supports the use of cellular oxygen and microcirculation on the surface of the skin. Helps drain lymph and swelling. At the same time, it balances and soothes the skin to achieve a stable effect.

EMS Tightening

EMS mode does not emit any light. The purpose of this mode is to send gentle electrical impulses to the muscles of the neck and décolleté, thereby strengthening them. This helps to firm the skin, improve skin tone and reduce wrinkles.

CLEANING AND MAINTENANCE

- **The device should be cleaned immediately after each use to prevent the accumulation of grease or other residues.**
- **Store the product in a dry place and at room temperature.**
- The product can be gently wiped with cold water, warm water, a soft cloth or a cotton swab.
- Do not clean the device directly under running water or immerse it, as this could damage the device.
- The product may be warm during use, but should never be hot.
- Do not use sharp, corrosive or volatile solvents, gas or other chemical cleaners that could damage the device.
- To maintain a long battery life, do not charge devices longer than necessary (3 hours).
- If the machine is not used for a long time, charge it at least once every 3 months and store it in the box.

RECOMMENDATION

When using the product, hold the silver bars on the sides of the device to fully realize the function of ion import and export. Please note that the product may not work properly if the strips are not held properly.

THE PRODUCT IS NOT INTENDED FOR:

Kids, women during pregnancy or menstruation, women with skin problems caused by hormonal imbalances, people with allergies, skin diseases, dermatitis, skin allergies and severe sunburns, as well as for people with skin damage, eczema and swelling. Also, the product must not be used by people with pacemakers.

The physically or mentally disabled are prohibited from using the product unless they are under the supervision and guidance of a guardian, as the guardian is responsible for the user's safety.

WARNING:

- **Do not use the device while charging the battery.**
- Before using the product, please read the user manual carefully and keep it properly so that other people who will use the product can also read the manual.
- If the charging cable is damaged, do not connect the product to the mains to avoid damaging other electrical devices. In case of damage or complaint, contact the manufacturer or the store where you purchased the product.
- If you experience pain, abnormal skin sensations or skin damage, stop using the product immediately and consult a doctor.
- Store out of reach of children. Prohibited to use by children under 6 years of age!
- Do not use the product near flammable materials (alcohol, gasoline, nail polish remover, etc.)
- Protect the product from moisture, do not put it in the dishwasher, washing machine, microwave oven or other household appliances.
- **Do not use the device in the shower or bath, or near any other water source** (water tanks, swimming pools, sinks, etc.).
- If the product is damaged by a careless drop or collision, do not use it again.
- Do not disassemble or repair the product at home without a professional. Contact the manufacturer or store where you purchased the product.
- Use only the supplied USB cable for power to ensure the required voltage is reached and product specifications are met.
- People with electronic medical devices, especially people with cardiac pacemakers, should only use the product after consulting a doctor.
- Do not use the device excessively long and often, and follow the recommended treatment time in the manual.

TECHNICAL PARAMETERS:

- **Product dimensions:** 50x92x135 mm
- **Battery voltage:** 3.7V
- **Frequency:** 5W, 24 thousand $\pm$ 1 thousand
- **Battery capacity:** 450 mAh
- **Charging time:** 2-3 hours

CHARGING:

Connect the USB cable to the product (at the bottom). Connect the USB cable to the computer or the adapter to the socket.

Images with descriptions are for illustration purposes only.

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SkyTech Corporation s.r.o., U libeňského pivovaru 2015/10,