Vögan[®]
Mammoth

Lady Three

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Vacuum skin cleaner with micro dermal abrasion

Vigan Mammoth Lady THREE digital

Dear customer, we would like to thank you for buying this device. Please read carefully this instruction manual before usage. Use this product only as described in this manual. Keep the instruction manual for later use. We recommend to keep the original packaging and receipt. If the device will be handed to another person, always attach the instruction manual.

General information


Vacuum microdermabrasion is one of the most popular aesthetic dermatology treatments. Thankful to this device, you can enjoy this procedure from your home. Lady THREE cleans and rejuvenates your skin and reduces facial wrinkles.

Microdermabrasion is a type of skin rejuvenation that prevents aging and reliably reduces facial wrinkles. This method effectively helps to remove blackheads (comedones), surface pigmentations, and facial skin irregularities after acne.

The special fine microscopic crystals of the attachment included in the package remove the micro-layers of dead cells from the epidermis. This treatment promotes cell renewal and regeneration of the skin surface. Its patency is significantly improved for other applied nourishing or curative remedies (nourishing creams). At the same time, the vacuum produced provides better blood circulation and the combination of all the effects of microdermabrasion gradually leads to increased formation of collagen and elastin fibers and thus to the much-desired rejuvenation. The result is a reinforcement of the basic layers of the skin and hence overall skin healing, which is smoother, tougher, more flexible and has a fresh look.

Microdermal abrasion does not cause any major contraindications, is painless and suitable for virtually all skin types. It perfectly meets the requirement of non-invasiveness and zero risk.

Vacuum lifting and peeling with Lady THREE works on the principle of vacuum. The skin is sucked into 4 different interchangeable inserts that are part of the device and due to the resulting vacuum, the slack skin is cleaned and stretched out. The fifth attachment is designed for light therapy.

5 attachments

USB cable

2 pcs of silicone
sealing rings

50 pcs of spare
micro-filters

Removable heads

LCD display

Control button

For switching of levels
and powering
the device down

First button push

– low suction intensity

Second button push

– medium suction intensity

Third button push

– strong suction intensity

Fourth button push

– blue light

Fifth button push

– red light

Sixth button push

– turning the device off

Device body

4 cm

16,5 cm

Time of use

Battery status

Operation indication

Suction intensity level
– low, medium and strong

Micro USB
charging connector

Vigan®
Mammoth

str 33

Microcrystalline peeling head

It contains large surface area microcrystalline particles that have absorbing properties, surface activity, etc. Without the risk of injury, it massages the skin and exfoliates, or flakes the dead skin. It helps clean and restore (stimulates the cell renewal process) the skin. It is not recommended to use this head around the eyes where the skin is too thin. Use only once a week if you have dry skin, twice a week if your skin is oily or combined. Always a maximum of 5 minutes.

Head with large round hole (larger suction intensity)

The larger suction force of this nozzle is intended primarily for vacuuming clogged pores (black dots), grease, dust, cosmetics residues, etc. Use the head carefully on thin skin. Use the head once a week for five minutes, more only depending on your skin condition. It is not recommended to use this head around the eyes, but it can be used for rougher and wrinkled skin parts, as well as for the body and elbows.

Silicone head with round hole (low suction intensity)

This head uses low suction power. Because silicone is relatively soft, the head can be used to absorb grease, dust, cosmetic residues, etc. It is designed for thin and delicate skin or normal skin.

Head with oval hole (low suction intensity)

The oval opening is designed for daily skin care. For example, for wrinkled corners of the mouth. After sucking in, quickly pull off the device and repeat the procedure. You can also use the face pulling technique from top to bottom, for example for faces, nose, mouth, etc.



Red and blue massage suction

Special design. Blue light - suitable for pore shrinking and soothing skin after sucking off dirt, acting on your skin antiseptic and anti-inflammatory. Red light - for rejuvenating effect and softer skin. Can be used on cheeks and nose. Let the red and blue light work on the skin at a distance of 3 cm for 3 to 5 minutes. Another option is to move straight on the skin, where a gentle massage occurs.

Micro filter interchange:

The key for the right treatment is to open the pores
the pores can be opened by hot shower, spa or steam

Main characteristic of Vigan Lady THREE

1. Device is designed for easy and effective skin care
2. 3 levels of suction intensity
3. 2 modes of light therapy (blue light) calms the skin and has antiseptic and anti-inflammation effect. (red light) - for rejuvenating effect and softer skin
4. Non-invasive, painless treatment for almost every type of skin
5. Vacuum lifting and peeling from the comfort of your home
6. Improves blood circulation and has anti-aging effect
7. Removes dead cells and supports its regeneration
8. Removes the impurities in the pores
9. Supports the microcirculation and strengthen the layers of your skin
10. Device stimulates the collagen and elastin production
11. Ergonomic shape for easy handling and manipulation
12. Easy to carry
13. Low weight and space saving design
14. USB charging
15. Charging time approx. 90 minutes
16. Up to 90 minutes of operation time

Safety instructions:

- Before use, make sure that the device and accessories do not show any visible damage. When in doubt, do not use it and contact your dealer or a professional service, to avoid injury.
- Dispose of packaging materials appropriately. Beware, for example, in the case of plastic bags, there is risk of suffocation. Make sure that children do not come into contact with packaging material (plastic bags, etc.).
- Keep and store the unit and the power unit in a dry environment.
- Before use, make sure that the microfilter is properly placed in the designated place, that it is dry so that no water or other liquids are absorbed by the device. At the same time, wipe water and other liquids out of the skin areas you want to treat, otherwise device will suck in the fluid, which blocks it and prevents it from running properly.
- Please remove any accessories (such as necklace, earrings, etc.) and contact lenses when using this device to prevent damage to the device and personal injury.
- Do not immerse the appliance in water! Do not use the device if your hands are wet!
- The device is intended for domestic use only. It is not possible, for example, to use it in beauty salons or by multiple persons at the same time
- Do not leave the device unattended while turned on
- Do not use or charge this device in a humid environment (for example, near a full bath tub, shower, pool, etc.)
- Never take the device into your hand if it falls into the water. Unplug the power cord immediately to avoid electric shock.
- Protect the device from shocks and be careful it will not to fall out of your hands.
- Protect the device from moisture, dust, chemicals, strong temperature fluctuations, and near heat sources (stoves, fireplaces, heaters).
- Be careful not to pinch or bend the power cord.
- Ensure that the power cord does not come into contact with water, steam, or other liquids.

- The device is not suitable for people under 18 years of age!
- Use the device exclusively on the human face (except eye parts, eyelids, ears, neck and mouth) of adults. Do not use the device on animals!
- The product is not intended for use on parts of the body where plastic operations have been performed, such as facelift or rhinoplasty (nose surgery). Furthermore, it is not recommended that the device is used by pregnant women and people suffering from heart disease.
- If you experience pain or other abnormalities on your skin or body, stop using it immediately.
- This device is not waterproof. It is strictly forbidden to vacuum with it water, immerse it in water or wash it directly under running water.
- **This device may be used by persons with reduced physical, sensory and mental abilities or lack of experience and knowledge when properly instructed and supervised in connection with the safe and proper use of this equipment. Cleaning and user maintenance should not be performed by children without supervision. Children should not play with the appliance.**
- Do not use accessories other than the original ones supplied by the manufacturer
- If the device has been damaged, for example, by dropping or dropping it, stop using it immediately.
- Do not use the device during charging.
- If the power cord is damaged, it may always be replaced with a cable of the same type to avoid possible injury and damage to the device.
- Do not use the appliance in water or in its immediate vicinity (bath, pool or shower).
- Do not use the device if the attachments are broken, cracked, dirty or otherwise damaged
- Avoid contact of any metal-containing material with the cable socket to avoid short circuits
- This device contains a magnet. Keep the device away from pacemakers.
- Persons with pacemakers or other electronic devices implanted in their bodies must not be in close proximity to the device. If a person is interested in using a magnet device or is in close contact with it, they must first consult a specialist doctor.
- Do not open the device or replace the rechargeable battery.

ATTENTION:

- Do not use the device for purposes other than those specified in this manual! Neither the supplier nor the seller can be held liable for damages resulting from improper use.
- Do not repair the device yourself. If damaged, contact your dealer or an authorized service center.

Principles to prevent skin irritation:

- Do not use Lady THREE more than twice a week
- Do not use an exfoliating microcrystalline attachment with other extensions in one day
- When using the device, the device must be in constant movement on the skin, never leaving the unit in one place for long periods of time, as this may cause irritation, injury and bruising.
- Do not use the device in conjunction with (peeling) creams or with mechanical and/or chemical peeling. This may result in skin injury and damage to the device.
- The device is intended for external use
- Occasionally, a reddening of the skin may occur, which is a normal condition that gradually disappears after the end of the process.

Do not use:

- For cracked skin, skin that is injured, reddened, irritated, bleeding, or if you have any other skin disease in the treated areas, do not use on rash, burns, blisters, ulcers, eczema, psoriasis, rosacea, cucerosis or spider veins, hypersensitive skin, inflamed sites, warts (could spread), active infection, or markedly active acne, if you have a big birthmark, tattoo, fire (nevus flammeus), permanent make-up in treated areas, or if the skin underwent a cosmetic procedure in the last two weeks, for example after chemical or mechanical peeling, laser or fractional laser treatment, or mesotherapy or other filler treatments, if you are taking skin-affecting medicines such as steroids or preparations for acne treatment.
- If you have scars in the area of treatment, it is also recommended that you consult your doctor or dermatologist before use.
- on the eye and mouth
- If you have skin cancer
- Immediately after sunbathing or sunburn
- If you are taking blood thinners
- If you have any of the following: multiple sclerosis, diabetes or autoimmune disease.

IF YOU HAVE ANY SKIN DISEASE, CONSULT A SPECIALIST BEFORE USING THE APPLIANCE.

IF YOU FEEL UNCOMFORTABLE WITH YOUR SKIN OR HAVE PROBLEMS WITH YOUR SKIN, CONSULT YOUR DOCTOR.

THE USE OF THE DEVICE DOES NOT REPLACE MEDICAL CONSULTATION AND TREATMENT. ASK YOUR DOCTOR BEFORE USING THE DEVICE.

Recharging

We recommend that you fully charge the device after each use.

- The appliance is intended for use, charging and storage at temperatures between 5°C and 35°C.
- To maintain battery life, fully charge the device at least every three months.
- If you have not used the unit for a long time or the battery has drained, recharge it before using it again.
- The unit cannot be used during charging.
- You can use your device even when the battery is not fully charged. In this case, however, the device's performance and endurance may not be optimal.
- Charging takes about 90 minutes and the device then runs up to 90 minutes depending on the intensity and type of attachments used.
- Make sure the device is switched off before charging.
- Plug the micro USB cable into the device socket and then into the USB input.

- When charging, the battery symbol appears on the display to indicate the charging process. You can recall this symbol at any time by pressing the control button on the display.
- When fully charged, the display shows, after pressing the control button, a static symbol of the charged battery, no longer signaling the charging process.
- When the battery is completely discharged, the device automatically turns off and the display stops to shine.

Before use

Always wash your face thoroughly and dry before the use of Lady THREE. Make sure that you removed the make-up, cream etc. We recommend steaming your face before the treatment- it will open your pores and improve the effect.

Use of device

Warning: 

Unpack the product from the box. Dispose of packaging materials appropriately. Beware, for example, in the case of plastic bags there is a risk of choking of children. Make sure that children do not come into contact with the packaging material.

Areas of face treatment

Nose: Move the device slowly from the bridge of your nose down to its tip and sides.

Temples: Move the device from outward to the hair line. Be careful to keep the device at least 1 cm from your eyes!

Areas of face treatment

Forehead: Move the device slowly from the center of the forehead towards the hair line. Make sure that you don't treat any spot more than once!

Cheeks: Move the device slowly from the sides of your nose towards your jaws.

Surrounding of the mouth:

Move the device from the center under the nose/chin to both left and right sides. Make sure that the device will never come in contact with your lips!

First use of the device

Before using the device for the first time, it is recommended to test the attachment on the inside of the forearm.

1. Attach the selected attachment to the top of the device as intended. With silicone sealing rings, a little force is needed and carefully mount with one-sided circular motion.
2. Turn the device by button push.
3. Press the button again to select the desired intensity from three levels. The fourth level is the activation of the blue light in conjunction with the skin calming attachment and the fifth level is the activation of the red light for rejuvenating effect and softer skin.
4. Before moving on to the face, move one stroke over the inside of the forearm and then wait about half an hour. If the eventual skin's reaction to this time subsides, you can go to your face.

Use - Tips:

Exfoliation of the dead skin on the face

1. Attach the microcrystalline attachment to the device. Wash and dry face.
2. Start from the lower jaw and move the device from the bottom up.
3. Place the device on your face and move it from the inside out.
4. Place the device on your forehead and move it from the center to the sides.
5. Place the unit on a T-shaped part of the face and move it from bottom to top.

Cleaning pores, grease, cosmetic residues, etc.

1. Attach a large round hole adapter or silicone ring adapter (depending on your skin type).
2. Point the nozzle to the skin parts you want to care for. Do not move the nozzle back and forth, but only in one direction. Do not use this mode for too long.

Care for fine lines

1. Attach the oval attachment to the device.
2. Use a device to treat the nose and lips. Adjust the suction force according to your needs. Take care of the part and pull it away quickly. Repeat the procedure several times for a maximum of 5 minutes.

Skin care with red and blue light

You can start this care right after cleaning. Attach the blue and red light attachment and move it very slowly. Do not use it in one place for too long to avoid skin redness and discomfort. Let the red and blue light act on the skin at a distance of 3 cm for 3 to 5 minutes. Another option is to go straight on the skin, where a gentle massage occurs.


First facial use of the device

For the first use on your face, we strongly recommend to start with the lowest level of suction so your skin can get used to this treatment. If there will be no difficulties after first application, you can increase the suction power.

1. Attach the device head that you wish to use on your face. Make sure that there are the silicone seal rings under the attachment. You can use circular movement to install the attachment properly.
2. Turn on the device
3. Press the button again to select the desired intensity from three levels. The fourth level is the activation of the blue light in conjunction with the skin calming attachment and the fifth level is the activation of the red light for rejuvenating effect and softer skin.
4. Place the attachment on the face skin so you can feel the suction. The attachment must be in permanent and full skin contact. It is not necessary to push the attachment too hard but only lightly. You should feel the attachment stick to your skin.
5. Swipe the device straight from the nose to temple and then wait about half an hour. If the eventual skin reaction has resolved within half an hour, you can safely use the device all over your face.
6. It is recommended to use the device twice a week and to keep the break between uses min. 48 hours.

After use

When you have finished using the device, apply for example a face mask or moisturizing cream, either day or night.

Follow the following precautions for at least 3 days after use:

- Never use any mechanical or chemical peeling.
- Postpone any planned electrolysis, wax treatment or depilation procedure in the face.
- Protect your skin from sun rays.
- You can apply makeup immediately after skin treatment.

Cleaning and maintenance

CAUTION: **Never immerse the device in the water!**

Never use rough wipes, abrasive scrubbers or aggressive chemical solution to clean the appliance! (petrol, toluene or acetone)

Cleaning:

1. Make sure that the device is turned off
2. Carefully remove the attachment
3. Clean the attachment under hot running water using dishwashing detergents. You can also use a soft brush to remove greasy dirt. In case of exfoliation, we recommend that you soak attachment for 10 minutes in the dishwashing solution and then clean with a soft brush.
4. Dry the cleaned attachment with a soft cloth so that no lint is left on it and put it back in the box, or put it back on the device.
5. Remember to replace the attachment filter if it is dirty.
6. Clean the device with a damp cloth and then dry with dry cloth.

Problem solving

Here are some of the most common problems and answers that may occur with your device.

The most common problems and questions – their solutions

This device is suitable for my skin?

Lady THREE is suitable for every skin type, for both oily, dry and mixed skin. It is also suitable for all skin colors. See the instructions below for more details and different warnings, or consult your doctor.

How long does it take to treat the skin with the Lady THREE

The treatment time for different parts of the face is different, but the complete facial treatment takes about 20 minutes.

Is Lady THREE technology safe

The device operates on the basis of safe microdermal abrasion technique. It only removes excess skin particles, enhances skin cell renewal, revives and regenerates skin surface. Its patency is significantly improved. The base layers of the skin are strengthened and hence the overall skin regeneration, which is smoother, tougher, more flexible and fresh in appearance.

How often to use the Lady THREE?

It is recommended to use the device twice a week provided that there is a break between uses for at least 48 hours. For highly sensitive skin, we recommend using the device max. once a week. Do not use an exfoliating microcrystalline attachment with other extensions in one day.

The device does not work,

Charge the device according to the section “Charging” in this manual. Make sure that there is no power failure and that there is a voltage in a power socket or device connected to the USB power cable. Check that the device is charging by pressing the button on the display. Make sure you press the power button correctly. If the device still does not work, contact a service center or importer.

Is the application of Lady THREE painful?

No, the application of the device should be painless. In very rare cases, mild burning and redness may occur in very sensitive skin. If you feel pain during application, we recommend lowering the intensity to a lower level.

Movement on skin is not smooth / weak or no suction / application result is worse than previous treatment?

Check microfilter for clogging. Replace or clean the filter. Check the attachment if it is clean. Check the peeling attachment if it is clean, or follow the “Cleaning and Maintenance” section. Make sure you feel the suction after removing the attachment. Turn on the device. When you place your finger on the opening, you should feel the suction. Failure to do so may indicate a malfunction. Please contact a specialist technician.

The device cannot be charged.

Check the connection of the device and whether the voltage is indicated on the device, or verify that the USB power cable is working. If the battery symbol still does not light up on the display and the device still does not charge, contact a service center or importer.

My skin reacts to the use of the device?

Your skin may react and become slightly reddened immediately after use. You can feel both warm and gentle tingling during use. These manifestations usually disappear within a few minutes. It may also cause your skin to be dry and full of flakes, and skin redness may take longer, or your skin may even swell. The reason may be a disproportionate pressure on the device or an increased treatment activity in the same area.

Do not apply pressure to the appliance and make sure that you do not treat the same area several times in succession or stay on the same place for a long time. You can also prevent these symptoms by applying a protective factor after using the device. The device has been tested by the manufacturer and has been confirmed to be safe to use as described in this manual.

Maintaining constant suction on the forehead is problematic.

Move the device more slowly in the forehead area to find the best contact with skin. Move the device slower in the forehead area.

Can I go to direct sun after treatment?

Protect your skin from direct sunlight for the first three days after treatment.

Can I perform other facial treatments after treatment?

Do not perform any mechanical, chemical peeling, electrolysis or facial epilation and do not use any depilatory agents for min. 3 days.

Technical specifications

Battery capacity: 500 mAh

Charging voltage: 5 V == 1 A

Operational voltage: 3,7 V

Rated power: 2,2 W

Protection class: III

Charging time: about 90 minutes

Operational time: up to 90 minutes

Color: white-grey

Weight netto: 154 g

Weight brutto: 432 g

Product dimensions netto: 165 x 40 mm

Dimensions brutto: 234 x 146 x 60 mm

Included in package

1x Gift box with padding for storing of the device and accessories

1x Unit Lady THREE

1x Microcrystalline peeling head

1x Head with large round hole

1x Silicone head with round hole

1x Head with oval hole

1x Red and blue massage suction

1x USB cable

2x Seal silicone ring

50x microfilter

1x User manual (CZ/ SK/ UK/ DE/ HU)

Disposal of old electrical and electronic equipment

This symbol on the product, its accessories, packaging or accompanying documents indicates that the product should not be treated as household waste. Please take this product to the respective collection point where will be carried out specialist recycling of electronic equipment. In the European Union and other European

countries there are collection points for used electrical and electronic equipment. By ensuring proper product disposal, you prevent potential negative impact for the environment and human health, which could result from improper handling of this product. Recycling of these materials also helps to preserve natural resources. For this reason, please do not throw the old electric and electronic equipment together with household waste. For more information please contact your local authority, any service company liquidating household waste or the store where you purchased the product. Details also refer to the relevant laws about waste, in the Czech Republic no. 185/2001 Coll. as amended. Furthermore, the websites www.elektrowin.cz, www.remasytem.cz and www.asekol.cz (disposal of electrical waste). The improper disposal of this waste can be in accordance with applicable laws, regulations and fines.

Warranty

We offer warranty for 24 months starting from the day of purchase by you. This warranty covers defect caused by defective material or error in workmanship. It does not cover defects caused by mechanical damage, natural disaster or other inevitable events, connecting to the wrong power grid voltage, modifications or service work by non-qualified or unauthorized personal, misuse or use in contradiction to this manual.

This device is equipped with long-life Li-ion battery. The 24 months warranty doesn't include the decreasing battery capacity caused either by the use or age. The capacity decreases depending on how often and how do you use/charge the battery.

All documents about your purchase and service repairs on any product keep for at least three years to ensure quality customer service.

Please remember, that this product is only for home use, not professional use.

If you find out, that this product does not suit you, it is not a reason for a complaint.

In the case of complaint please submit receipt.

In the case of malfunction of your product, please contact the dealer, you bought it from, or service at the following address:

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