GALVANIC FACE IRON 5 IN 1

USER MANUAL

PROFESSIONAL TREATMENT AT HOME

**DEAR CUSTOMER,
THANK YOU FOR PURCHASING THIS PRODUCT.**

PLEASE READ THIS USER MANUAL CAREFULLY AND KEEP IT FOR FUTURE REFERENCE. THIS PRINTED MANUAL CONTAINS THE USER MANUAL IN 4 LANGUAGES. FOR OTHER LANGUAGE VERSIONS, PLEASE SCAN THE QR CODE.

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EN GALVANIC FACE IRON 5 IN 1

PALSAR7 - ABOUT THE BRAND

DEAR CUSTOMER, THANK YOU FOR PURCHASING THIS PRODUCT. PLEASE READ THIS MANUAL CAREFULLY BEFORE USE AND RETAIN IT FOR FUTURE REFERENCE.

PALSAR7 IS A BRAND SPECIALIZING IN INNOVATIVE PRODUCTS FOR PROFESSIONAL SKIN CARE FROM THE COMFORT OF YOUR HOME. IT OFFERS MODERN TECHNOLOGIES THAT SUPPORT SKIN HEALTH AND BEAUTY WITHOUT THE NEED TO VISIT A BEAUTY SALON. THE BRAND EMPHASIZES QUALITY, SAFETY AND EASE OF USE OF ITS PRODUCTS. THANKS TO THIS, YOU CAN INDULGE IN PROFESSIONAL CARE ANYTIME, ANYWHERE.

ABOUT THE PRODUCT

PALSAR7 5-in-1 Facial Galvanic Device combines several modern technologies in one elegant device, allowing you to enjoy advanced skincare from the comfort of your home.

It utilizes **galvanic technology, micro vibrations, heating, and LED light therapy**. These functions help support a **cleaner, fresher, and more radiant** appearance of the skin. They also promote better **absorption** of **serums and creams**, making your daily skincare routine more effective.

Thanks to several specialized modes, you can easily **adapt the treatment** to your skin's **current** needs and provide it with the care it deserves.

SIDE STRIPS

Keep your **fingers on both side silver strips** during use.

EFFECTS

5 skincare technologies with **ions**, **red** and **blue LED light**, **heat**, and **micro vibrations**.

SUITABLE WITH COSMETICS

Use with skincare products such as **creams** and **serums**.

EFFECTIVE

It helps promote **smoother**, **firmer**, and more **radiant** skin.

MODES

1. WARM CLEAN MODE

The Warm Cleanse Mode uses gentle heat, vibrations, and ionic technology **to more thoroughly remove impurities** from the skin's surface and pores. The heat helps open the pores and prepare the skin for more effective cleansing. At the same time, **the ions help** draw impurities and makeup residue onto the makeup remover pad. After treatment, the skin appears **cleaner, fresher, softer, and naturally radiant**.

- **Frequency:** Use 2–3 times a week.
- **Instructions:** Place a makeup remover pad with makeup remover on the head and turn on the device. While using the ionization function, make sure to keep your skin well-moisturized at all times. Place your fingers on both side silver strips and gently cleanse the skin using circular motions. After finishing, remove impurities from the face with a clean cotton pad.

2. PHOTOTHERAPY MODE / RED LED LIGHT: 640 ± 10 NM

The Phototherapy Mode is designed to promote **hydration, nourishment, and healthy-looking skin**. When used in combination with a toner, serum, or cream, it helps these **skincare products absorb better**. It promotes **skin elasticity** and helps **reduce the appearance of fine lines**. After treatment, the skin appears **smoother, more nourished, and better hydrated**.

- **Frequency:** Use in the morning and evening. If you have sensitive skin, reduce the number of treatments.
- **Instructions:** Apply toner, cream, or essence to cleansed skin and turn on the device. Place your fingers on both silver strips on the sides and gently glide the device across your face in an upward motion.

3. SKIN LIFTING MODE BLUE LED LIGHT: 450 NM ± 10 NM

The Skin-lifting Mode helps promote **more comfortable and intensive skin care**. A pleasant warmth helps promote **microcirculation** in the skin, promote **skin elasticity**, and contribute to a firmer appearance. After treatment, the skin appears **smoother, refreshed, and revitalized**.

- **Frequency:** You can use it daily.
- **Instructions:** Apply your favorite face mask to your face and turn on the device. Place your fingers on both silver strips on the sides and gently glide the device across your face in an upward motion.

MODES

4. EYE CARE MODE


The Eye Care Mode is designed for **gentle treatment of the sensitive area** around the eyes and helps eye cream or essence absorb better into the skin. It supports the skin's **natural hydration and nourishment**, leaving it soft, fresh, and with a healthy glow. The intensity of this mode is constant.

- **Frequency of treatment:** You can use it daily, in the morning or evening.
- **Instructions:** Apply your favorite skincare products, such as eye cream or serum, to your face and turn on the device. Place your fingers on the two silver strips on either side and use the device to gently massage the area around your eyes. If needed, you can use the intensity control button to turn on the red light, which helps reduce the appearance of fine lines.

INSTRUCTIONS FOR USE

Before first use, charge the device for at least 2 hours.

1. **Clean the facial device** before each use.
2. Use only on **clean skin**. For better glide and effectiveness, apply a suitable **cosmetic product** according to the selected mode, such as a makeup remover, toner, cream, essence, or eye cream.
3. Press and hold the lower button  for **2 seconds** to turn the device **on or off**. Then press the button once to switch between functions.
4. After selecting a mode, you can adjust the vibration intensity in selected modes using the **upper** button on the handle. The selected mode and intensity are displayed on the device screen.
5. During use, it is necessary **to touch both silver conductive strips located on the sides of the handle**.
6. Move the device **slowly and gently** without applying excessive pressure. On the face, guide the device using **circular motions**. Be especially gentle around the eye area.
7. The device does not switch off automatically. Please turn it off manually after use.

MAINTENANCE AND CLEANING

- Before cleaning, make sure that the device is **switched off and disconnected** from the power supply.
- Clean the device before each use according to the instructions below.
- The device should **be cleaned immediately** after each use to prevent the buildup of oil, dirt, or other residues.
- Store the product **in a dry place at room temperature**.
- To clean the device, use a **soft cloth** or a **cotton pad** lightly moistened with **lukewarm or cold water**.
- **Do not clean the device under running water** and do not immerse it in water, as this may damage the device.
- The product may become warm during use, but it should never become hot.
- **Do not use** sharp objects, corrosive or volatile solvents, gas, or other chemical cleaners that could damage the device.
- To help maintain battery life, do not charge the device longer than necessary (**3 hours**).
- If the device is not used for an extended period, charge it at least once every **3 months** and store it in its box.

TECHNICAL SPECIFICATIONS:

- Rated Voltage: **3.7 V / 0.6 A**
- Battery Capacity: **350 mAh**
- Charging Time: **2-3 hours**
- Red LED Light: **640 ± 10 nm**
- Blue LED Light: **450 ± 10 nm**

NOTICE

This product is not intended for:

- **children, women** during **pregnancy**, or women during menstruation
- individuals with **skin problems** caused by hormonal imbalance
- individuals with **allergies, skin conditions, dermatitis**, or **severe sunburn**
- individuals with **damaged skin, eczema**, or **swelling**
- people with a **pacemaker** or other **electronic medical implants**
- individuals with reduced physical, sensory, or mental capabilities may use the device only **under the supervision** of a person responsible for their safety

CHARGING:

Connect the USB cable to the port on the side of the device. Connect the other end of the USB cable to a computer or to a **standard wall power adapter**.

PACKAGE CONTENTS:

- 1× PALSAR7 Massage and cleaning facial iron 5 in 1
- 1× USB Cable
- 1× User Manual
- 1× Plastic head attachment

DIMENSIONS:

- 168 x 43 x 50 mm

WARNING:

- Do not use the device **while the battery is charging**.
- Before using the product, carefully read the user manual and keep it in a safe place so that **other users can also read it if necessary**.
- If the charging cable is damaged, do not connect the product to a power source to avoid damaging electrical equipment. In the event of damage or a warranty claim, contact the manufacturer or the retailer where the product was purchased.
- If you **experience pain, unusual skin sensations, or skin damage, stop using the product immediately** and consult a physician.
- **Keep out of reach of children.**
- If you wish to use the device near a **dental prosthesis**, exercise increased caution.
- Do not use the product near flammable materials (alcohol, gasoline, nail polish remover, etc.).
- **Protect the product from moisture. Do not place it** in a dishwasher, washing machine, microwave oven, or other household appliances.
- Do not use the device in a **shower or bathtub, or near any water source** (water tanks, swimming pools, sinks, washbasins, etc.).
- If the product is **damaged by a fall or impact, do not use it again**.
- **Do not disassemble or repair** the product yourself. Contact the manufacturer or the retailer where the product was purchased.
- For power supply, **use only the USB cable included in the package**, ensuring the required voltage and compliance with product specifications.
- Do not use the device excessively long or frequently, and follow the recommended treatment times stated in the manual.

POSSIBLE OPERATIONAL ISSUES

The device cannot be turned on:

1. Try pressing and holding the power button for a few seconds.
2. If the battery is discharged, recharge the device.

The device does not vibrate:

1. Make sure you are touching both silver strips on the handle.
2. Use a moistened cotton pad, serum, toner, or another suitable cosmetic product. The device may not function properly on dry skin.

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Wear and tear of the product and its components caused by their common use (discharged batteries, loss of the original color of the product, decreased capacity of accumulator, wear and tear of mechanical or displaying components, wear and tear of supply cables, etc.) shall be subject neither to defect warranty nor quality warranty in the moment of take-over.

Furthermore, complaints about product parameters, not included in the operating instructions or in any commercial-technical documentation of Seller, related to the product, shall not be subject to quality warranty in the moment of take-over.

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