medivon

User manual

Cosy Vibra

PL | EN | DE | FR | ES | PT | IT | CZ | SK
SL | EE | LT | LV | HU | RO | RU | NB | FI
EL | HR | SV | DA | NL | TR | BG | AR | UK

PL

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
EN

1. Setup and operation

1.Secure the cushion in place by wrapping the back elastic straps over the car seat or chair.
2.Connect the adapter cable to the corresponding cable at the side of the massager.

3. Plug the home adapter into an outlet or the car adapter into your car's power receptacle.
4. Turn on the device using the controller.
5. When finished place the controller in the pouch at the side of the seat massager and unplug the adapter from the outlet or receptacle.


• Do not setup or use the massager in a bathroom or similary wet/damp areas as such an environment may reduce the lifespan of the massager mechanism.
• Ensure the massager is not exposed to direct sunlight or placed near sources of heat.

Note: This product produces a warming heat. It is not intended to be used as a therapeutic heating pad. Please allow a few minutes to reach maximum temperature. The cushion reaches a maximum temperature of 51°C (124°F) after approx. 7 minutes.

2. Safety Instructions

Please read the following instructions carefully before putting your massage device into use to ensure trouble-free operation and optimum efficiency. Please retain these operating instructions for further use!

• Do not wet, do not use pins, never remove cover.
• This items is NOT A TOY. Close supervision is necessary when this appliance is used by, on, or near children or disabled persons.
• This appliance should never be left unattended when plugged in.
• Any possible repairs may only be carried out by authorized expert staff.
• Improper use and unauthorized repairs are not permitted for safety reasons and lead to loss of warranty.
• Never touch the power plug with wet hands.
• Please avoid contact of the device with water, high temperatures and direct sunlight.
• Do not use any damaged cables, plugs or loose sockets.
• Never operate under blankets where the air opening may be blocked.
• If plugs or cords are damaged, they must be replaced by the manufacturer, a service representative or by qualified personnel.
• In case of malfunction, disconnect immediately from the mains.
• Do not use if you have skin disorders, open wounds,or swollen or inflamed areas.
• Misuse or incorrect use excludes any liability for damage.
• Do not use while you sleep.
• In order to avoid excessively stimulating the muscles and nerves, the recommended massage time should not be exceeded 60 minutes at a time.
• Each massage – even a hand massage – must be refrained from during pregnancy or if one or more of the following complaints are present in the massage area: recent injuries, thrombotic diseases, all kinds of in flammations and swellings, and cancer. It is recommended to consult a doctor before having a massage for the treatment of ailments and illnesses.
• If you are dependent on electrical aids e.g. pacemakers, please be sure to consult your doctor for medical advice before taking a massage.
• The packaging material supplied should not be used as a toy.

3. Controller instructions

1. Body zone; indicator lights;

2. Intensity indicator lights (LO, MED, HI);

3. On/Off;

4. Intensity selector;

5. Push for Neck/Upper back massage;

6. Push for Mid-back massage;

7. Push for Lower Back/Lumbar massgae;

8. Push for Leg / Thigh massage;

9. Mode indicator lights (M1 - M5);

10. Push to turn heat on / off;

11. Mode selector - choose from 5 pre-programmed massage modes.

Massage
Each massage area may be selected separately. Select high, medium or low intensity. For use at home.

Heat
Press the "HEAT" button to turn the heat function on or off. A built-in thermostat controls the temperature within a safe range. The heat may be used with or without massage.

NOTE: This product produces a warming heat. It is not intended to be used as a therapeutic heating pad. Please allow a few minutes to reach maximum temperature. Mat reaches maximum temperature of 60°C after approx. 7 minutes.

Built-in Timer
The cushion has a built-in timer that deactivates the cushion after 15 minutes. After 15 minutes of use please wait 15 minutes and simply turn cushion on again.

4. Care and Cleaning

- Disconnect from power source before cleaning.
- Do NOT place in any liquid to clean. Clean fabric cover by wiping with a damp sponge, using a mild detergent or suds. Pat dry with clean, dry cloth.
- Do NOT use abrasive cleaners, brushes, gasoline, kerosene, glass or furniture polish, paint thinner or water hotter than 49°C to clean the unit.

5. Storage

Place the unit in its box, or in a safe, dry, cool place. Avoid contact with sharp edges or pointed objects that might cut or puncture the fabric. Do not wrap power cord around the unit. Do not hang the unit by the power cord.

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Used electrical and electronic equipment and batteries and accumulators must not be mixed with other waste. This is indicated by the crossed-out bin symbol on the product packaging. Electrical and electronic equipment and batteries and accumulators may contain hazardous substances, mixtures and components that are harmful to the environment and human health. Improper handling or damage to them may result in damage to health or environmental pollution. Proper handling of used electrical and electronic equipment and batteries and accumulators in the household contributes positively to the reuse and recovery, including recycling, of used equipment. Used electrical and electronic equipment and batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and has performed the same functions as the equipment sold. A distributor, when delivering household equipment to a purchaser, is obliged to take back used household equipment free of charge at the place of delivery of such equipment, provided that the used equipment is of the same type and has performed the same functions as the equipment delivered. A distributor running a retail outlet with a sales area within the meaning of art. 2 item 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m^2 devoted to the sale of household equipment, is obliged to accept, free of charge, in that outlet or in its immediate vicinity, used household equipment, none of the external dimensions of which exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market after 13 August 2005.

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Short Declaration of Conformity
The manufacturer declares under its sole responsibility that the Cosy Vibra product complies with Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, and Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU as regards the list of restricted substances.
The full text of the Declaration of Conformity is available at the following website: www.medivon.pl/deklaracje

Importer: Armare S.A.
ul Grzybowska 87,
00-844 Warszawa, Poland
Made in China
Marka: Medivon
E-mail: serwis@medivon.pl
www.medivon.pl
SN:01265908275917113

Producent: XIAMEN COMFIER TECHNOLOGY CO.,LTD
2-5/F, No.65 Sunban South Road, Jimei District, Xiamen, China.
E-mail: sales@comfier.com

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Shortened warranty terms
This Guarantee is provided by Armare S.A., a company governed by Polish law, with its registered office in Warsaw, ul. Grzybowska 87, 00-844, Poland. You can find us in the official Polish register of companies under the KRS number: 0000730590, NIP: 5213833614 and REGON: 380721860. We guarantee that the equipment you have purchased from us (hereinafter referred to as the "Product") will work properly. If something goes wrong for a reason inherent in the Product, this warranty protects you. The warranty starts from the day you receive the Product and lasts: 24 months, if you are a consumer or a natural person running a business and the purchase of the Product is not of a professional nature for you. Otherwise, the warranty period is 12 months. If the law in your area requires a longer warranty, the warranty period will be extended to meet local law requirements.

In addition, if the Product does not comply with the contract, you have the right to legal protection from and at the expense of the seller. This warranty does not affect these rights - it does not limit, suspend or exclude them.

Detailed warranty conditions, how to submit complaints, what happens when you find a defect, what are exceptions to the warranty, when the warranty expires or is extended, and what are the rules regarding the processing of personal data, you can find on our website: http://medivon.pl/warranty